

# BIO BM200

## Portable Mini Vein Finder

### Brief Introduction:

Vein finder BM200 is easy to detect veins at the back of the hands through transmission of LED red light, and it can assist medical staff to observe sub-cutaneous veins of hands, and significantly improve the success rate of the injection.

It is suitable for hospitals, clinics, homes as well as any other places to find the veins for injection.

### Product Features:

- small in size, light in weight, easy to carry
- adjustable light, wide application range
- sensitive switch, safe and power saving
- fit human body engineering, more comfortable grip
- rechargeable lithium battery
- Automatic light sensor for turn/off(optional)



### Method of Operation:

- turn on the rotary switch
- hold the light bulb with palm. Now the vein finder is sending light.
- rotate the rotary switch, adjust the strength of light, the veins are visible (more dark than other tissues).
- after vein puncture, turn off the rotary switch.



Wuxi Biomedical Technology Co., Ltd.

[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu

Tel: 86-510-85572789 Fax: 86-510-85572789

Email: sales@cnbiomedical.com